

package deals horizontals

Cleaver Scientific now offers a comprehensive range of Combination Packages to address any budget or horizontal electrophoresis application.


Mini Horizontal with Power Supply



Ordering Information

| | |
|-----------------------------|---|
| MSMINIDUO-MINI150 | Mini Horizontal Electrophoresis Package, includes: - |
| | <ul style="list-style-type: none"> • MSMINIDUO - horizontal unit with 7x7cm and 7x10cm gel trays and combs (Page 6) • MINI150V - 150V, 600mA, 60W mini power supply (Page 39) |
| MSMINIDUO-MINI150KIT | Mini Horizontal Electrophoresis Package with electrophoresis starter kit, includes: - |
| | <ul style="list-style-type: none"> • MSMINIDUO-MINI150 • CSL-LOADDYE - 10x Bromophenol Blue Loading Dye, 1ml (Page 15) • CSL-MDNA-1kb - 1kb DNA ladder (Page 15) • CSL-AG100 - 100g agarose (Page 15) • CSL-TBEP - Powdered Tris-Borate-EDTA Running Buffer - 10 pouches (Page 15) |

Midi Horizontal with Power Supply



Ordering Information

| | |
|-----------------------------|---|
| MSMIDIDUO-MINI300 | Midi Horizontal Electrophoresis Package, includes: - |
| | <ul style="list-style-type: none"> • MSMIDIDUO - horizontal unit with 10x7cm and 10x10cm gel trays and combs (Page 7) • MINI300V - 300V, 400mA, 60W mini power supply (Page 39) |
| MSMIDIDUO-MINI300KIT | Midi Horizontal Electrophoresis Package with electrophoresis starter kit, includes: - |
| | <ul style="list-style-type: none"> • MSMIDIDUO-MINI300 • CSL-LOADDYE - 10x Bromophenol Blue Loading Dye, 1ml (Page 15) • CSL-MDNA-1kb - 1kb DNA ladder (Page 15) • CSL-AG100 - 100g agarose (Page 15) • CSL-TBEP - Powdered Tris-Borate-EDTA Running Buffer - 10 pouches (Page 15) |

Choice Horizontal with Power Supply



Ordering Information

| | |
|------------------------------|--|
| MSCHOICETRIO-CS300 | Choice Horizontal Electrophoresis Package, includes: - |
| | <ul style="list-style-type: none"> • MSCHOICETRIO - horizontal unit with 15x7cm, 15x10cm and 15x15cm gel trays and combs (Page 8) • CS-300V - 300V, 700mA, 150W midi power supply (Page 39) |
| MSCHOICETRIO-CS300KIT | Choice Horizontal Electrophoresis Package with electrophoresis starter kit, includes: - |
| | <ul style="list-style-type: none"> • MSCHOICETRIO-CS300 • CSL-LOADDYE - 10x Bromophenol Blue Loading Dye, 1ml (Page 15) • CSL-MDNA-1kb - 1kb DNA ladder (Page 15) • CSL-AG100 - 100g agarose (Page 15) • CSL-TBEP - Powdered Tris-Borate-EDTA Running Buffer - 10 pouches (Page 15) |

Maxi Horizontal with Power Supply



Ordering Information

| | |
|------------------------|--|
| MSMAXIDUO-CS300 | Maxi Horizontal Electrophoresis Package: - |
| | <ul style="list-style-type: none"> • MSMAXIDUO - horizontal unit with 20x10cm and 20x20cm gel trays and combs (Page 9) • CS-300V - 300V, 700mA, 150W midi power supply (Page 39) |